



Trusted by Hospitals. Loved by Everyone.

# Manufacturer Information for Use

## SONO Disinfecting Wipes

### 20 / 50 / 80 COUNT SOFTPACK - 80 / 100 COUNT CANISTER

---

**Use SONO Disinfecting Wipes in:** Ambulances, Hospitals, Medical Offices, Nursing homes Health Clubs, Schools, Airplanes, Emergency Vehicles.

**Use SONO Disinfecting Wipes on washable hard, nonporous surfaces of:** Countertops, Desks, Walls, Computer keyboards.

**Non-Critical Medical Equipment Surfaces:** Examining Tables, Stethoscopes, Ultrasound transducers and equipment.

**Other hard nonporous surfaces made of:** Glass, Metal, Plastic, Stainless steel, Upholstery, vinyl and plastic.

**To Open Package:**

Slowly peel back the label. Remove wipe as needed. To retain moisture, reseal by firmly pressing seal back on the pouch.

**Directions for Use:**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Do not use as a diaper wipe or for personal cleansing. This is not a baby wipe! Do not flush in toilet. For plastic and painted surfaces, spot test on an inconspicuous area before use. This product will not harm most surfaces: including acrylic, sealed fiberglass and vinyl.

**To Clean and Disinfect Virucidal\* Directions:**

Wipe hard, non-porous surface with wipe until surface is visibly wet. Keep surfaces wet for 4 minutes. Allow surface to remain wet for 4 minutes. For Norwalk Virus, Feline Calicivirus, Norovirus and Rotavirus, treated surfaces must remain wet for 10 minutes. Let air dry. If surfaces are visibly dirty, clean first with another wipe before disinfecting. Discard used wipe in trash. Do not flush.

**Mold and mildew stain control,** SONO Ultrasound Wipes will effectively inhibit the growth of mold and mildew and the odors caused by them when applied to hard, nonporous surfaces. Wipe surface and let it remain visibly wet for 4 minutes. Repeat treatment every seven days, or more often if new growth appears.

**To Clean and Deodorize:** Wipe surface with wipe. Let air dry. Discard used wipe in trash. Do not flush in toilet.

**\*KILLS HIV AND HBV AND HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS**

**PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 [HIV-1] [associated with AIDS] or Hepatitis B Virus [HBV] or Hepatitis C Virus [HCV].

---

---

## **Special Instructions for Cleaning and Decontamination against HIV-1 or HBV or HCV on Surfaces/Objects Soiled With Blood/Body Fluids:**

**Personal Protection:** Clean-up must always be done wearing protective gloves, gowns, masks and eye protection.

**Contact time:** Leave surface wet for 1 minute for HIV-1. Use appropriate contact time as listed on label for all other viruses and bacteria claimed.

**Cleaning Procedure:** Blood and other body fluids containing HIV or HBV or HCV must be thoroughly cleaned from surfaces and objects before application of SONO Ultrasound Wipes. This product can be used for this purpose.

**Disposal of infectious material:** Blood, body fluids, cleaning materials and clothing must be disposed of according to local regulations for infectious waste disposal.

This product is for use on hard nonporous surfaces not conducive to treatment by immersion or excess liquid against Avian Influenza A. This product is not to be used as a terminal sterilant / high level disinfectant on any surface or instrument that [1] is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or [2] contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**SANITIZING AND DEODORIZING:** To sanitize and deodorize hard, non-porous, non-food contact surfaces; wipe surface, using enough wipes to allow treated surface to remain wet for 15 seconds for *Staphylococcus aureus* [Staph] [ATCC 6538], *Klebsiella pneumoniae* [ATCC 4352], and *Campylobacter jejuni* [ATCC 29428]. Let surface dry. For visibly soiled surfaces, clean excess dirt first.

### **STORAGE AND DISPOSAL**

Store in original container in areas inaccessible to children. Keep securely closed. Nonrefillable container. Do not reuse or refill this container. Discard empty container in trash or recycle if available.

### **SONO Disinfecting Wipes kills:**

#### **Bacteria:**

*Pseudomonas aeruginosa* [Pseudomonas], *Salmonella enterica* [Salmonella], *Staphylococcus aureus* [Staph], *Acinetobacter baumannii*, *Burkholderia cepacia*, *Campylobacter jejuni* [Campylobacter], *Corynebacterium ammoniagenes*, *Enterobacter aerogenes*, *Enterococcus faecalis* [Enterococcus], *Enterococcus faecalis*- Vancomycin Resistant [VRE], *Escherichia coli* [E. coli], *Escherichia coli* O157:H7, ESBL *Escherichia coli* – [Extended spectrum beta-lactamase producing E. coli], *Klebsiella pneumoniae* [Klebsiella], *Klebsiella pneumoniae* (NDM-1) [New Delhi metallo-beta-lactamase 1], *Legionella pneumophila*, *Listeria monocytogenes*, *Salmonella schottmuelleri* [Salmonella], *Salmonella typhi* [Salmonella], *Serratia marcescens* [Serratia], *Shigella dysenteriae* [Shigella], *Staphylococcus aureus* - Community Associated Methicillin-Resistant [CA-MRSA] [NRS384] [USA300], *Staphylococcus aureus* - Community Associated Methicillin-Resistant [CA-MRSA] [NRS123] [USA400], *Staphylococcus aureus* - Methicillin-Resistant [MRSA], *Staphylococcus aureus* - Multi-Drug Resistant [resistant to tetracycline [Tc], penicillin [Pc], streptomycin [Sm] and erythromycin [Em] and susceptible to chloramphenicol [Cm] in vitro], *Staphylococcus aureus*-Vancomycin Intermediate Resistant – [VISA], *Streptococcus pyogenes* [Strep] [a cause of scarlet fever], *Vibrio cholerae*

---

---

**Viruses:**

\*SARS-Related Coronavirus 2 [SARS-CoV-2] [cause of COVID-19], \*Hepatitis B Virus [HBV], \*Hepatitis C Virus [HCV], \*Herpes Simplex Virus Type 1, \*Herpes Simplex Virus Type 2, \*HIV-1 [AIDS Virus] , \*Human Coronavirus, \*Influenza Type A Virus / Brazil [Influenza], \*Influenza A Virus H1N1, \*2013 Influenza A Virus (H7N9), \*Norwalk Virus [Feline Calicivirus] [Norovirus], \*Respiratory Syncytial Virus [RSV], \*Rotavirus, \*SARS Associated Coronavirus [SARS] [cause of Severe Acute Respiratory Syndrome], \*Vaccinia Virus

**Animal Viruses:**

Avian Influenza Virus H3N2, Avian Influenza Virus H5N1, Canine Distemper Virus, Feline Calicivirus , Newcastle's Disease Virus, Pseudorabies Virus

**PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

**Caution.** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

**FIRST AID**

**IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling poison control center or doctor or going for treatment.



Trusted by Hospitals. Loved by Everyone.