

SONO[®] Disinfecting Wipes

TRUSTED BY HOSPITALS.
LOVED BY EVERYONE.



Hospital Grade Disinfectant Wipes

- Included on the EPA (U.S. Environmental Protection Agency) "List N" for disinfectants for Coronavirus (COVID-19).
- Destroys pathogens including MRSA, HIV and staph; inhibits the growth of mold and mildew as well as the odors caused by them.
- Non-abrasive, non-damaging formula has been tested compatible with all major brands of ultrasound equipment.
- Effective on hard, washable, non-porous surfaces in many healthcare settings, as well as at home, in the gym or office, or when traveling.
- Top-quality fabric with antibacterial, antiviral and antifungal formula; quaternary ammonium based.
- Textured to pick up particles and quickly remove dirt, spots and smudges while killing 99.9% of germs, bacteria, and viruses.



Distributed by:
NEXT Medical Products Company
customerservice@nextmedicalproducts.com

SONO® Disinfecting Wipes

TRUSTED BY HOSPITALS. LOVED BY EVERYONE.



Effective against Coronavirus (COVID-19), MRSA, HIV, and staph.
Kills 47 different pathogens. Contact us for a full listing.

FOR ULTRASOUND

Most general disinfectant wipes can cause damage to ultrasound equipment, leading to costly repairs or replacement. SONO Wipes have been tested to clean, disinfect, deodorize and sanitize ultrasound equipment -- including transducers and LCD/touch screens -- without causing harm or damage. Protect your investment with SONO Wipes.

- ✓ Tested on major brands of ultrasound equipment.
- ✓ Kills 47 pathogens including MRSA, HIV and Staph.
- ✓ Safe to handle without gloves.
- ✓ Hospital grade disinfectant that won't damage equipment.
- ✓ Helps avoid costly damage to equipment such as delamination of transducer heads and yellowing and cracking of cables.
- ✓ Free of phosphates, ammonia, solvent and bleach; low alcohol level.
- ✓ Convenient - no need for multiple cleaning products. Can be used on non-porous surfaces as well as transducers, cables, lenses, LCD/ touch screens and system consoles.

Using the wrong wipe can cause serious damage.

Approximate cost to hospital: \$4,500
Touch screen panels became unresponsive after cleaning with hydrogen peroxide wipes

Delamination caused gel to seep inside.

Specialty Transducers can cost upwards of \$25,000
Alcohol and solvent can dry and crack cables, exposing wires

Replacement Cost \$9,000-\$15,000
Damage can lead to patient shock and image quality issues

FOR OTHER APPLICATIONS



In addition to ultrasound equipment, SONO Wipes can be used to disinfect EKG machines, stress units, patient monitoring units, mammography plates and other medical devices. SONO Wipes are also ideal for other washable, hard, non-porous surfaces in a multitude of settings such as at home, in the gym, in the office and when traveling.

For more information or to place an order, call 1-800-458-4254
or email: customerservice@nextmedicalproducts.com

Packaging Information:

Soft Pack
P/N 4018

Distributed by:



Wipe size: 7" x 10"
Packaged 50 per pack,
12 packs/case



Canister
P/N 4032

Wipe size: 7" x 8"
Packaged 80 per canister,
6 canisters/case

NEXT Medical Products Company
45 Columbia Road, Branchburg, NJ 08876 USA
nextmedicalproducts.com

Made in the USA
Meets OSHA's Bloodborne
Pathogen Standard 1910.1030

Patent Pending
SONO is a registered trademark of
Advanced Ultrasound Solutions, Inc.