

Distributed by



Combi[®] Laryngoscopes

All-in-one disposable laryngoscope set
with permanently attached blade and
handle, plus battery



No Compromise
Professional Instruments



Metal Max™ Disposable Combi® Laryngoscope Sets

Distributed by:



The Combi is an all-in-one disposable laryngoscope set with a **permanently attached** metal blade and handle, plus 6 volt battery, packaged in a color-coded, Test-in-Bag protective sheath.

- Blade:** Strong, medical grade stainless steel blade that extends through the hook-on for maximum strength
- Handle:** Robust, lightweight full metal aluminum
- Battery:** 6 volt with long battery life -- over 5 hours of continuous illumination

- The ready to use Combi eliminates compatibility issues and helps reduce SKUs.
- Saves time and labor related to managing and cleaning separate blades and handles.
- Eliminates the need to open blades and handles to insure that you have a pair that works.
- Supplied in unique, color-coded, Test-in-Bag packaging to allow immediate confirmation that the light works prior to opening the package.
- Battery can be removed quickly for recycling at end of life.
- Three year battery shelf life from date of manufacture.
- Available in a complete range of sizes including a lever tip configuration.



Part No.	Size/Description	Blade Length mm	Blade Length inches
LMMCB90MAC0	MAC 0	61.00	2.40
LMMCB90MAC1	MAC 1	75.00	2.95
LMMCB90MAC2	MAC 2	93.00	3.66
LMMCB90MAC2A	MAC 2 with Adult handle	93.00	3.66
LMMCB90MAC3	MAC 3	110.00	4.33
LMMCB90MAC3S	MAC 3 with Stubby handle	110.00	4.33
LMMCB100LH3	MAC 3 Left Handed	110.00	4.33
LMMCB90MAC35	MAC 3.5	122.00	4.80
LMMCB90MAC4	MAC 4	135.00	5.31
LMMCB90MAC4S	MAC 4 with Stubby handle	135.00	5.31
LMMCB100LH4	MAC 4 Left Handed	135.00	5.31
LMMCB90MAC5	MAC 5	156.00	6.14
LMMCB90MIL00	MIL 00	43.50	1.71
LMMCB90MIL0	MIL 0	55.00	2.16
LMMCB90MIL1	MIL 1	80.50	3.17
LMMCB90MIL15	MIL 1.5	100.00	3.94
LMMCB90MIL2	MIL 2	133.00	5.24
LMMCB90MIL2A	MIL 2 with Adult handle	133.00	5.24
LMMCB90MIL2S	MIL 2 with Stubby handle	133.00	5.24
LMMCB90MIL3	MIL 3	170.00	6.72
LMMCB90MIL4	MIL 4	182.00	7.17
LMMCB90LT3	Lever Tip 3	110.00	4.33
LMMCB90LT4	Lever Tip 4	135.00	5.31
LMMCB90CP00	Cardiff Pro CP 00	55.00	2.17
LMMCB90CP0	Cardiff Pro CP 0	71.00	2.81
LMMCB90CP1	Cardiff Pro CP 1	77.00	3.03

Combi sizes 00, 0, 1 and 2 are equipped with a pediatric handle unless noted



- PROACT INNOVATIONS**
- Introduced world's first fully metal disposable blade and handle
 - Introduced first fixed set with permanently attached blade and handle
 - First manufacturer to cover the light pipe to prevent back light



Manufactured by:



PROACT Medical Ltd
9-13 Oakley Hay Lodge, Great Folds Road
Corby, Northants, NN18 9AS UK
www.proactmedical.co.uk



North and South American Distribution by:
NEXT Medical Products Company
45 Columbia Road, Branchburg, NJ 08876 USA
www.nextmedicalproducts.com
800-458-4254, 908-722-4549, Fax 908-218-3715
customerservice@nextmedicalproducts.com