

ULTRASOUND GEL PACKET WARMER



The industry's first gel packet warmer, carefully designed for both clinical satisfaction and patient comfort.

- ~ Choice of three temperature settings
- ~ Digital display with real time temperature sensing
- ~ Holds 35-40 ultrasound gel packets
- ~ Tinted lid allows inside viewing
- ~ Easy to use, user-friendly interface
- ~ Proprietary warming controls maintain a consistent temperature throughout the day
- ~ Built-in safety feature helps prevent temperatures from reaching above the selected setting
- ~ Easy to wipe clean and disinfect

Tabletop Gel Packet Warmer

for Clear Image Singles®
Ultrasound Gel Packets

and other brands of gel in single packs

Settings:

Low 95°F (35°C)

Medium 100°F (38°C)

High 105°F (41°C)



Part Number:
07-01



Clear Image Singles is a registered trademark
of NEXT Medical Products Company, LLC

©2018 NEXT Medical Products Company

NEXT Medical Products Company
45 Columbia Road
Branchburg, NJ 08876 USA
www.nextmedicalproducts.com
800-458-4254, Fax 908-218-3715
customerservice@nextmedicalproducts.com

70207A