

# The Cardiff Pro Pediatric Laryngoscope Blade

Distributed by:



## *Finally! A Pediatric Standardization Option*

Proact Cardiff Pro Blades are specially designed\* (similar to Miller and Phillips blades) with a 5 degree downward angled blade (upward angled by 5 degrees when in use). It has a flat atraumatic tip that enables easier tongue control in pediatric patients.

Correctly formed blade shape and atraumatic tip where appropriate, no sharp edges.

Shielded optics provide properly focused output without distracting backlight.

Disposable, yet has a full-length spatula similar to a quality reusable blade.

Great for children's hospitals and pediatric departments. Available in three sizes -- 00, 0 and 1

\*The Cardiff Pro originated as a result of a meeting among pediatric anesthesiologists in Great Britain who joined together in Cardiff, Wales to discuss the optimal pediatric blade. Rather than having to choose from the more than two dozen blade options in the market, the Cardiff Pro was manufactured by Proact to combine the best features from all the other options -- thus allowing, standardization on ONE blade type.



GREEN SYSTEM



COMBI® LARYNGOSCOPE  
attached blade and handle set



CONVENTIONAL

Manufactured by:



PROACT Medical Ltd  
9-13 Oakley Hay Lodge, Great Folds Road  
Corby, Northants, NN18 9AS UK  
[www.proactmedical.co.uk](http://www.proactmedical.co.uk)



North and South American Distribution by:  
NEXT Medical Products Company  
45 Columbia Road, Branchburg, NJ 08876 USA  
[www.nextmedicalproducts.com](http://www.nextmedicalproducts.com)  
800-458-4254, 908-722-4549, Fax 908-218-3715  
[customerservice@nextmedicalproducts.com](mailto:customerservice@nextmedicalproducts.com)