

# Clear Image Singles® STERILE Ultrasound Gel Packets

*The original single patient gel packet  
is available in a sterile version.*

- **Gel is Free of Added Dyes and Perfumes.** Clear gel reduces the likelihood of skin sensitization and staining. It may also help prevent premature degradation of the transducer.

- **Gentle Formula.** Gel is water-based, water-soluble and contains no parabens, surfactants, oils, alcohol or silicone.

- **Control the Flow of Gel.** The built-in dispensing spout allows clinicians to direct the gel as it is removed from the packet -- from a dab to a glob.



- **Quality.** As specialists in ultrasound gel production, NEXT Medical is committed to producing the highest quality gel with the fewest bubbles possible. Made in the USA.

- **Convenient Packaging.** Each packet contains 20 grams of sterile ultrasound gel, overwrapped in a peel-apart outer pouch. Packaging Options:

- 04-048 48/box, 6 boxes/case
- 04-048D 48/dispensing carton, 6 cartons/case
- 04-048-300 Bulk Pack 300/case, no peel pouch



- **Easy to Use.** Peel apart the outer pouch



to remove the inner gel packet. Then, tear along the perforation. Gel can be dispensed easily through the built-in dispensing spout.

- **Performance.** Our viscous formula stays put when you need it to, yet spreads easily. It helps make secure contact between the transducer and the patient's body, preventing air pockets and increasing acoustic transmission. Designed for ultrasound applications where a sterile gel is desired.

**Try Clear Image Singles® STERILE  
Ultrasound Scanning Gel.**

**Call 1-800-458-4254 for more information.**

CE 2460 ISO 13485:2016



45 COLUMBIA ROAD, BRANCHBURG, NJ 08876 USA  
www.nextmedicalproducts.com  
Customer Service: 1-800-458-4254  
Email: customerservice@nextmedicalproducts.com