

ISO 7376 : 2020



## Metal Max™ COMBI™ 50 Laryngoscopes

All in one soft grip Combi™ system with pre-fitted alkaline batteries



6V Alkaline Battery

Latex Free

High strength Polymer handle with a soft grip moulding provides superior rotational and longitudinal grip.



Non-user detachable hinged surgical stainless steel blades with correct profile and form.

New **Metal Max Combi 50** offers premium quality

### Test in Bag Packaging Ready in Seconds

Our innovative and user friendly packaging allows no touch, full activation and function testing of the Laryngoscope in the bag.

More than just an illumination and switch test, this crucially allows you to check the full fit and function of ALL parts, including the LED and power source in their functional position.

Easy to identify colour coding, profile and large sizing information ensures you spend less time looking for the device you need when you need it.

Super bright, 25 Lumen LED light source, provides a minimum of 10,000 Lux\* across the range.



LED integrated into the blade

\*Light intensity measured at the tip of each blade using Lutron LX-1102 Light Meter and ambient shielding cylinder.

METAL MAX™

COMBI 50

One Piece Combi™ System



# Metal Max COMBI 50

ISO 7376:2020 Compliant

## Disposable One Piece System with Pre-Loaded 6V Batteries

- Pre fitted, non user detachable, Stainless steel blade with plastic handle. Supplied pre-loaded with 6V alkaline batteries - Ready for use in seconds.
- Unique, bayonet clip battery cap allows battery removal in less than 3 seconds
- NO TOUCH battery removal, allows clean and safe battery recycling / disposal
- Colour coded sachets with full functional 'test in bag' capability for peace of mind
- Ergonomic, moulded soft grip, provides excellent longitudinal, rotational grip and comfort
- 6V alkaline power source and flicker-free LED provide super bright illumination
- Full stainless steel blade design, providing maximum strength
- Atraumatic tip and fully formed blade profile with no sharp edges
- Full length ET Tube Guide surface
- Latex free
- 3 Year shelf life



Size 1 and smaller blades are supplied on a Paediatric Handle



Strength of stainless steel through the Hinge



Pre-installed Alkaline 6V Batteries

in a more cost effective system.

**NO TOUCH** battery removal, ready in...



◀ ...Seconds



Macintosh



Cardiff Pro



Miller



All the following are supplied in boxes of 10

Part No. Box of 10	Blade Design	Handle Design	Blade Size	ISO 7376 Intended Patient Use	Blade Length (mm)	Blade Width (mm)
LMMCB50MAC0	Macintosh	Paediatric	0	Neonate	61.0	8.5
LMMCB50MAC1	Macintosh	Paediatric	1	Small Child	75.0	11.5
LMMCB50MAC2	Macintosh	Adult	2	Child	93.0	12.5
LMMCB50MAC3	Macintosh	Adult	3	Adult	110.0	13.5
LMMCB50MAC4	Macintosh	Adult	4	Large Adult	135.0	14.0
LMMCB50MIL00	Miller	Paediatric	00	Premature Infant	44.0	10.5
LMMCB50MIL0	Miller	Paediatric	0	Neonate	55.0	10.5
LMMCB50MIL1	Miller	Paediatric	1	Small Child	81.0	10.5
LMMCB50MIL2	Miller	Adult	2	Child	131.0	13.0
LMMCB50MIL3	Miller	Adult	3	Adult	171.0	13.0
LMMCB50CP00	Cardiff Pro	Paediatric	00	Premature Infant	55.0	8.0
LMMCB50CP0	Cardiff Pro	Paediatric	0	Neonate	71.0	9.0
LMMCB50CP1	Cardiff Pro	Paediatric	1	Small Child	77.0	10.0

All devices manufactured by PROACT™ Medical are fully compliant with the relevant requirements of the European Medical Device Regulations MDR 2017/745.

All of PROACT's products are CE manufactured/FDA Approved to exceed the minimum requirements of any global standards. Our operations are supported by a Quality Management System certified to meet the criteria of ISO 13485:2016 and we expect the same from our suppliers.

We are proud of our environmental credentials and utilise a system which has been certified to the International standard ISO 14001. We constantly work to minimise our environmental impact, handling our business in a way that is committed to preserving, protecting and enhancing the environment.

TRADEMARKS: The following Trademarks belong to PROACT Medical Ltd: PROACT™, Metal Max™, COMBI™.



PROACT Medical Ltd  
9-13 Oakley Hay Lodge,  
Great Folds Road,  
Oakley Hay Business Park,  
Corby,  
Northants, NN18 9AS, UK

Tel: +44 (0)1536 461981  
Fax: +44 (0)1536 461837  
e-mail: sales@proactmedical.co.uk

proactmedical.co.uk

North and South American Distribution by:



NEXT Medical Products Company  
45 Columbia Road, Branchburg, NJ 08876, USA  
Tel: 908-722-4549 or 800-458-4254 Fax: 908-218-3715  
e: customerservice@nextmedicalproducts.com  
website: www.nextmedicalproducts.com

## Distributor Details

©PROACT Medical Ltd 2023. All rights reserved.  
PRO21DOC03-v2-2023-04

@PROACTMedicalLtd



@PROACTMedical



PROACT-Medical-Ltd-UK



@PROACTMedicalLtd



PROACT Medical

